上海科技大学
ShanghaiTech University

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
Loadlock放气

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
Loadlock抽真空

Vent Pick Wafer Place Wafer LL Init Start

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The screenshot displays the LEUVEN LIVELL-ICP control interface. At the top, there are status indicators for Warn, Alarm, and LL.Sts, along with buttons for Standby and LogOut. The main area is divided into two panels: 'ICP-RecipeList' on the left and 'ICP-StepList' on the right. The 'ICP-RecipeList' panel shows a table of recipes with columns for RecipeName, Creator, CreateDate, Modifier, ModifyDate, and Remark. A red dashed box highlights the 'RecipeName' and 'Remark' fields, with a label 'Recipe名称等信息'. Below this, a 'Recipe中包含的step' section shows a list of steps with columns for Index, StepName, and Remark. The 'ICP-StepList' panel shows a table of steps with columns for StepName, Creator, CreateDate, Modifier, ModifyDate, and Remark. A label 'Step列表' is placed below this table. At the bottom of the interface, there are buttons for 'New', 'Update', 'Load', 'Unload', and 'Delete'. The bottom navigation bar includes 'Process', 'Edit', 'Log', 'Users', 'Maintenance', 'Vent', 'Pick Wafer', 'Place Wafer', 'LL Init', and 'Start'. The bottom status bar shows 'LEUVEN-instruments Co. Ltd' and 'V2.0.0'.

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
The screenshot displays the control interface for the LIVELL-ICP system. At the top, there are status indicators for Warn (0) and Alarm (0), and a speaker icon. The PMSSts and LLSts status is shown as Standby. The main title is "LIVELL-ICP". On the right, there are buttons for "Init", "Abort", and "Shut Down" (highlighted with a red dashed box). Below these are fields for "User" (1) and "Level" (Maintenance), along with a "Logout" button. The central part of the interface features a schematic diagram of the system with various components labeled, including valves (YV1-1S, YV1-1F, TV1-1, QV1-1, LYV), pumps (DP1-1, LMP), and gauges (3E-003, 15.7 mTorr, 0.4 mTorr). On the left, there are control panels for various gases: CHF3 (2-100sccm), SF6 (3-150sccm), Ar (6-300sccm), O2 (6-300sccm), CF4 (6-300sccm), N2 (4-200sccm), Cl2/N2 (4-200sccm), and BCl3/N2 (4-200sccm). On the right, there are sections for "RF Top Information" and "RF Bottom Information", each with "On" and "Off" buttons and various power and setpoint parameters. Below these are controls for "Chiller" (OutPort T(°C) 40.0, Set T(°C) 40.0), "Cooling He" (0-100Torr, 1.8 Torr, 0.0 sccm), and "Electrode" (HeVA, HeVB, HeVC, HeVD). At the bottom, there are buttons for "Vent", "Pick Wafer", "Place Wafer", and "LL Init". The bottom left corner shows "Edit", "Log", "Users", and "Maintenance" buttons. The bottom right corner shows "TP1-1" and "DP1-1" speed and temperature controls.

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