

### Extraction from adherent cell lines

- (i) Change the medium of the cell plate(s) 2 h before metabolite extraction.
- (ii) Aspirate the medium completely.
- (iii) Put the plates on dry ice and add 2 ml of 80% (vol/vol) methanol (cooled to  $-80^{\circ}\text{C}$ ).
- (iv) Incubate the plates at  $-80^{\circ}\text{C}$  for 20 min.
- (v) Scrape the plates on dry ice with cell scraper.
- (vi) Transfer the cell lysate/methanol mixture to a 2-ml conical tube on dry ice **and Sonicate 10 min.**
- (vii) **After incubate 1h at  $-20^{\circ}\text{C}$ ,** centrifuge the tube at  $14,000g$  for 5 min at  $4-8^{\circ}\text{C}$  to pellet the cell debris.
- (viii) Transfer the metabolite-containing supernatant to a new 2-ml conical tube on dry ice.
- (ix) Add  $500\ \mu\text{l}$  80% (vol/vol) methanol ( $-80^{\circ}\text{C}$ ) to the pellet in a 2-ml tube and vortex for 1 min .
- (x) Spin the tubes at  $14,000g$  for 5 min at  $4-8^{\circ}\text{C}$ .
- (xi) Transfer the supernatant to a 2-ml conical tube on dry ice.
- (xiii) SpeedVac/lyophilize to a pellet using no heat.

Dried metabolite samples can be stored at  $-80^{\circ}\text{C}$  for several weeks.