

组织样品

Normalized by weight

Make sure <10% error of weight across different samples before extraction

1. Cut tissue on dry ice and record the tissue weight, make sure <10% error across different samples; (干冰上切割称重或液氮中研磨后称重)
2. 20 mg tissue + 200 μ L H_2O to homogenized, homogenize for three cycles (each cycle: 5500rpm for 20s, repeat three times) with cooled N_2 gas flow from liquid N_2
3. 200 μ L 匀浆液 + 800 μ L MeOH/ACN (v:v, 1:1)
4. Vortex 30 s, sonicate 10 min (4 $^{\circ}C$ water bath)
5. Incubate 1 h at -20 $^{\circ}C$ (fascinate protein precipitation)
6. Centrifuge 15 min at 13000 rpm and 4 $^{\circ}C$
7. Take supernatant and evaporate to dryness at 4 $^{\circ}C$ using a vacuum concentrator
8. Reconstitution: 100 μ L of ACN / H_2O (v:v, 1:1)
9. Vortex 30 s, sonicate 10 min (4 $^{\circ}C$ water bath)
10. Centrifuge 15 min at 13000 rpm and 4 $^{\circ}C$, placed in -80 $^{\circ}C$