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
Company Name : LUVA SYSTEM      Temperature Compensation : NO  
Model Number : Auto Four-point Probe 280 SI      SPC Option : YES  
Software Version : 5.3R6  
Current User : su      Program Description : Select Measurement Type

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Enter the process definition in the boxes provided from Page 1 to Page 6 :  
Edit Process : 12      Measurement Type : Circular pattern test

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Lower Data Plotting Limit :       Ohm/sq

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Upper Data Plotting Limit :       Ohm/sq

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Preview  
Print  
Set Skip List  
Set Leak Sites

OK      Cancel

	Company Name : LUYA SYSTEM	Temperature Compensation : NO
	Model Number : Auto Four-point Probe 280 SI	SPC Option : YES
	Software Version : 5.3R6	
	Current User : su	Program Description : Set Process Parameters

Enter the process definition in the boxes provided from Page 1 to Page 6 :  
Edit Process : 12 Measurement Type : Circular pattern test

Page 1	Page 2	Page 3	Page 4	Page 5	Page 6
--------	--------	--------	--------	--------	--------

OPTIONAL	Operator Name :	
OPTIONAL	Fab Name :	
OPTIONAL	Lot Identification :	
OPTIONAL	Equipment Information :	
OPTIONAL	Device Number :	
FIXED	Probe Head Serial No. :	
FIXED	Probe Head Type :	
FIXED	填写Probe探头类型: A/B	Probe Head Selection : Not Available
FIXED	填写探头校准因子, 或者从探头校准程序自动写入	Probe Head Correction Factor : 1
FIXED	Probe Pin Spacing :	1 mm

Preview  
Print  
Set Skip List  
Set Leak Sites

OK Cancel

	Company Name : LUYA SYSTEM	Temperature Compensation : NO
	Model Number : Auto Four-point Probe 280 SI	SPC Option : YES
	Software Version : 5.3R6	
	Current User : su	Program Description : Set Process Parameters

Enter the process definition in the boxes provided from Page 1 to Page 6 :  
Edit Process : 12 Measurement Type : Circular pattern test

Page 1	Page 2	Page 3	Page 4	Page 5	Page 6
--------	--------	--------	--------	--------	--------

FIXED	Wafer Size :	(for rectangular wafer, type W x H) 150 mm
	硅片尺寸100 125 150 200, 特殊尺寸请按照WxH输入	
OMIT	Flat / Notch Wafer :	NOTCH
	硅片形状	
FIXED	Max Test Diameter :	(for rectangular wafer, type W x H) 140 mm
	最大测量直径, 需比Wafer size至少小6mm	
FIXED	Mask Circle Diameter :	0 mm
FIXED	Number of Test Points :	(for rectangular wafer, type W x H) 65 points
FIXED	Gain Control :	AUTO
FIXED	Geometric Correction :	OFF (Enter Correction Factor above)

Preview  
Print  
Set Skip List  
Set Leak Sites

OK Cancel

	Company Name : LUYA SYSTEM	Temperature Compensation : NO
	Model Number : Auto Four-point Probe 280 SI	SPC Option : YES
	Software Version : 5.3R6	
	Current User : su	Program Description : Set Process Parameters

Enter the process definition in the boxes provided from Page 1 to Page 6 :  
Edit Process : 12 Measurement Type : Circular pattern test

Page 1	Page 2	Page 3	Page 4	Page 5	Page 6
OMIT	Edge Compensation :		OFF		
FIXED	Temperature Compensation :		OFF		
FIXED	Data Pre-sort Limit :		<small>(0-15, 0=disable)</small> 0	%	
					数据初筛选, 偏离中心值多少的去掉
FIXED	Data Sorting Limit :		<small>(0-6, 0=disable)</small> 3	Sigma	
					数据筛选正态分布区间
FIXED	Save Data on Disk :		YES		
					自动保存测试数据
FIXED	Autodraw Contour Map :		YES		
					每次测试完, 是否自动画图

Preview  
Print  
Set Skip List  
Set Leak Sites

OK Cancel

	Company Name : LUYA SYSTEM	Temperature Compensation : NO
	Model Number : Auto Four-point Probe 280 SI	SPC Option : YES
	Software Version : 5.3R6	
	Current User : su	Program Description : Set Process Parameters

Enter the process definition in the boxes provided from Page 1 to Page 6 :  
Edit Process : 12 Measurement Type : Circular pattern test

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Preview  
Print  
Set Skip List  
Set Leak Sites

OK Cancel

Company Name : LUYA SYSTEM      Temperature Compensation : NO  
 Model Number : Auto Four-point Probe 280 SI      SPC Option : YES  
 Software Version : 5.3R6  
 Current User : su      Program Description : Set Process Parameters

Enter the process definition in the boxes provided from Page 1 to Page 6 :  
 Edit Process : 12      Measurement Type : Circular pattern test

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FIXED    Auto SECS-II Data Transfer (serial) : OFF

OMIT    Notch Reference Pattern :

FIXED    Set Max Current Range : (-1:high to 4:low)

OMIT    Measure Leakage Current : NO

FIXED    Auto Data File Save : NO  
 自动导出文件存储，选No

Preview  
 Print  
 Set Skip List  
 Set Leak Sites

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


Please enter Wafer identification below : [Press F2 for name same as previous]

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


First Configuration Measurement

Second Configuration Measurement

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(Ohms/sq)	(Best Gain)	(Ohms/sq)	(Best Gain)
0 to .175	LOW	175 TO 700	HIGH
.175 TO .70	HIGH	700 TO 1750	LOW
.70 TO 1.75	LOW	1750 TO 7000	HIGH
1.75 TO 7.0	HIGH	7000 TO 17.5K	LOW
7.0 TO 17.5	LOW	17.5K TO 70K	HIGH
17.5 TO 70	HIGH	70K TO 175K	LOW
70 TO 175	LOW	175K TO 800K	HIGH

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Type	Application	Material	Tip Radius @ Contact Area Radius	Force
Type-A	Metal films, ITO, amorphous Si, poly silicon, SOI, bulk substrates	Tungsten Carbide	25 um or 45um	90-200g
Type-B	General Purpose: medium to thick epitaxial layers, diffused layers, metal films, implants 2-10 microns	Tungsten Carbide	100um	90-200g
Type-M	epitaxial, diffused and ion implanted layers (500 A to 5 um), especially very shallow ion implanted layers	Tungsten Carbide	300um @ 120um	90-200g
Type-N	Very shallow implants	Tungsten Carbide	500um	90-200g
Type-C	III-V compound semiconductors, e.g. GaAs, InAs	Osmium	100um	60-150g
Mercury	Ultra shallow implants, ultra thin layers	Mercury		~0
Custom	Customized tips	Various	Various	

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